

Digital Dry Bath

Description

Welso Digital Dry Baths WDDB series offer digital control over both temperature and time, providing accurate and uniform temperature distribution on the thermal blocks. It is molecular biology, histology, clinical, environmental and industrial laboratories. Our digital dry bath are available in three configurations with capacity for up to four blocks.

- Real time and temperature display
- Over- temperature protection provides more safety
- Compact, fits on benchtop and saves space
- High quality aliminum heating blocks



Specification

Model	WDDB100	WDDB200	WDDB400
Blocks	1	2	4
Temperature	RT+5°C~150°C		
Heating time(20°C-150°C)	≤ 30 m		
Temperature Stability at 40°C	°C		
Temperature	°C		
Temperature Stability at 120°C	± 1 °C		
Temperature Display Accuracy	0.1°C		
Time Setting	0-99h59m		
Max. Temperature	150°C		
Power	200W	400W	600W
Dimension(mm)	260×220×100	260×220×100	360×220×100

Block Option

Order #	Block		
		DDBB010	8*26mm
DDBB001	48*6mm	DDBB011	6*28mm
DDBB002	48*7mm	DDBB012	3*40mm
DDBB003	24*10mm	DDBB013	48*0.2ml
DDBB004	24*12mm	DDBB014	48*0.5ml
DDBB005	24*13mm	DDBB015	24*1.5ml
DDBB006	12*15mm	DDBB016	24*2.0ml
DDBB007	12*16mm	DDBB017	
DDBB008	12*19mm	DDBB018	96*0.2ml
DDBB009	12*20mm	DDBB019	96 well Elisa Plate
			Customized Blocks upon to request